



Pressure and  
Temperature Measurement

## New adapter system for sterile processes

Klingenberg, February 2012.

**WIKA has now developed a process adapter system for measuring instruments with hygienic connections, which meets the high requirements of sterile engineering processes. It consists of a pressure gauge or transmitter with a mounted adapter and is certified to 3-A Sanitary Standards.**

The flexible, modular system enables the connection of a wide range of aseptic process connections such as clamp, threaded, VARIVENT® or NEUMO®. All parts are manufactured from 316L stainless steel. The O-ring for sealing the process (optional), available in EPDM or FKM, is approved by FDA, USP Class VI and 3-A 18-03.

With this adapter system, customers can switch to measuring instruments with different process connections in a very short time. Therefore they can reduce long downtimes and inventory costs.

Number of characters: 858  
Key words: Process adapter

**Manufacturer:**

WIKA Alexander Wiegand SE & Co. KG  
Alexander-Wiegand-Straße  
63911 Klingenberg  
Tel +49 9372 132 – 0  
Fax +49 9372 · 132 – 406  
E-Mail Vertrieb@wika.de  
[www.wika.com](http://www.wika.com)

press release

**WIKAI company photograph:**

Fig. top: Pressure transmitter model S-11 with G1 hygienic connection

Fig. centre: Process adapter VARIVENT®

Fig. bottom: Clamp process adapter



**Edited by:**

WIKAI Alexander Wiegand SE & Co. KG

André Habel Nunes

Marketing Services

Alexander-Wiegand-Straße 30

63911 Klingenberg/Germany

Tel. (+49) 9372 132-8010

Fax (+49) 9372 132-8008010

E-Mail [a.habel-nunes@wika.de](mailto:a.habel-nunes@wika.de)

[www.wika.de](http://www.wika.de)